



Zimmer®
Drill
Stop Kit



A smarter solution for controlling drilling depth.



Zimmer® Drill Stop Kit. An open and shut case for an easier, more cost-efficient drill stop solution.

Beauty and simplicity, together in one cool package.



Drill stop technology simplified.

With securely fitted, length-specific stops, the *Zimmer Drill Stop Kit* is designed to simplify the drilling process and save chair time for you and your patients.

Time savings.

- Eliminates the need to carefully follow drill depth markings and for frequent pauses to examine osteotomy depth

Ease of use.

- “Pick & go” drill stop application delivers easy engagement onto the drill
- Stop diameter color-coding correlates to *Zimmer Instrument Kit System* color-coding for easy identification
- Drill length labeling aids in easy identification
- Unique Multi-tool allows for easy pick-up of stops from the kit, and stop removal from the drill
- Two measurement guides aid in drilling depth and stop location verification for increased safety

Cost-efficient design.


- Reusable, sterilizable titanium material eliminates the need for continuous investment in replacement stops



Drill stops for use with 17/22mmL drill

Columns are color-coded by implant diameter and organized by drill diameter

Rows are organized by implant lengths




Drill stops for use with 11/16mmL drill

Multi-Tool
Multi-tool aids in pick-up and removal of stops


Sample Sequence

Osteotomy for a 3.7mmD x 13mmL Tapered Screw-Vent® Implant, using a 17/22mmL Driva™ Drill.



SV2.3DN
2.3mmD
Drill

Step 1: From the 13mmL implant row, select the stop for a 2.3mmD Pilot Drill.



SV2.8DN
2.8mmD
Drill


Step 2: From the same row, select the stop for a 2.8mmD drill (final for soft bone) or skip to the stop for a 3.4/2.8mmD Step Drill (final for dense bone).



TSV3DN
3.4/2.8mmD
Drill

Taking depth control to new heights.

In any osteotomy, depth control and patient safety are major concerns. You want a drill stop system that is easy to understand and easy to use, and most important, delivers what you need to perform implant surgery with confidence. That's where the new *Zimmer Drill Stop Kit* comes in. With convenient "pick & go" simplicity, easy-to-follow color-coding and organization, and a host of other time-saving features, the *Zimmer Drill Stop Kit* is designed to rise above all other depth control solutions.



Step 1:

Select the appropriate drill stop.

Step 2:

Verify drilling depth by using the Drill Depth Guide.

Step 3:

Create the osteotomy to a predetermined depth.

Step 4:

When done, simply disengage the drill stop using the Multi-tool or by hand.

Step 5:

Before placing the drill stop back into the kit, verify the stop's location in the kit by using the Drill Stop Guide to ensure proper placement.

For more information about our Products, Regenerative Materials and Educational Opportunities, contact us:

1900 Aston Avenue
Carlsbad, CA 92008-7308, USA

In the U.S. 1 (800) 854-7019

To fax an order 1 (888) 225-2483

Outside the U.S. +1 (760) 929-4300

Australia +61 (0)2 9950 5434 or 1 (800) 241 916

Canada +1 (905) 567-2073 or 1 (800) 265-0968

France +33 (0)1 45 12 35 35

Germany +49 (0)761 1 56 47 0

Israel +972 (0)3 6124242

Italy +39 043 855 5573

Spain +34 93 846 05 43

Nederland:

Implacom BV

Paleisweg 5 - 3886 LC Garderen

Postbus 57 - 3886 ZH Garderen

T. 0577 461927

E. info@implacom.nl

I. www.implacom.nl

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsletterSignUp.aspx



www.zimmerdental.com

