


AccuraMesh™

Product Information

CustomGraft Solutions™ for Guided Bone Regeneration



AccuraMesh Products

CustomGraft Solutions For Guided Bone Regeneration

Bone regeneration procedures and techniques have advanced the way that hard tissue is rebuilt prior to implant placement. Bone blocks and bone graft particulates are available from a choice of human, animal, or synthetic origin. However, most of the available products are pre-shaped and need to be adapted to the defect site. Their surfaces require manipulation in order to achieve ingrowth of blood vessels and vascularization of the graft material.¹

Now, with the widespread use of digital technologies in modern dentistry, clinicians have the advantage of using 3D imaging systems, 3D-printed guides, or patient-specific prosthetic components to achieve customized solutions. Pre-designed bone grafting materials can be manufactured based on a CBCT/CT scan of the patient's defect area.

Introducing Zimmer Biomet AccuraMesh! Available in either Titanium or PEEK, Zimmer Biomet AccuraMesh is designed using a fully digital workflow. Data from 3D medical imaging devices combined with modern Computer-Aided Design (CAD) software and state-of-the-art Computer-Aided Manufacturing (CAM) processes result in high-quality customized medical devices for guided bone regeneration procedures.²



Your Benefits

Zimmer Biomet AccuraMeshes Have The Following Features:

1

Two Material Options Available

- PEEK
- Titanium

2

Pre-Planned Screw Positions

- Reliable fixation
- Reducing the risk of touching sensitive anatomical structures³

3

Sterile Packaged – Ready For Use

- 10⁻⁶ Sterility Assurance level^{4,5}

4


Space Maintaining

- Protects and secures bone graft particulates for undisturbed healing⁶

5

Reduced Surgery Time And Morbidity

- Additional manual adjustment of the defect and of the customized meshes is seldom required, allowing for reduced surgery time and reduced morbidity^{7,8}



Two Different Raw Materials, Two Different Products

Titanium AccuraMesh

Made of Surgical Grade Titanium Alloy by Selective Laser Melting Technology^{2,4}

- Proven biocompatibility⁹
- Strong enough to maintain space⁴
- High precision

Electropolished Surface²


- Reliable mesh removal¹⁰
- Reduced bacterial adhesion¹¹
- Enhanced soft-tissue cell adhesion¹²

Pre-Planned Hooks

- Membrane fixation

Pre-Planned Screw Positions

- Reliable fixation



AccuraMesh Products

Two Different Raw Materials, Two Different Products


PEEK AccuraMesh

Made of Implantable Grade Polyether Ether Ketone Filament (PEEK) by Fused Deposition Modeling Technology

- Proven biocompatibility¹³
- Designed for long-term implantable medical devices⁵
- Residual metal ions according to ISO 10993-18 less than 0.5 ppm¹³

Smooth Surface

- Reduced bacteria adhesion¹⁴



Pre-Planned Screw Positions

- Reliable fixation

Please Note The Following Information

Imaging

Patient Preparation

- Remove temporaries and metal restorations, where possible
- Position patient in stable position

Imaging Requirements

- In general, most CT/CBCT devices are suitable
- Recommended slice thickness: 0.2 to 0.75 mm
- Gantry angle: 0°
- Open bite scan
- Please ensure that high-contrast imaging is achieved

Scan Data

- Do not use data compression
- Data must be provided in DICOM format only*
- Transfer the files using Zimmer Biomet Dental Upload website: www.zbd-fileupload.com

Planning And Design

Design Draft

- You will receive by e-mail (I) a 3D-PDF file of the designed customized mesh and the defect site and (II) design & validation form
- To open the PDF files, Adobe Acrobat Reader is required
- Adjustments can be made at any time before final approval

Surgery

- Select proper flap design and soft-tissue management to ensure tension-free soft-tissue closure¹⁵⁻¹⁸
- Despite precise planning, the products may not fit as expected and minor manual adjustments may be required

Ordering Information

Item #	Description
TICMS	Titanium AccuraMesh Standard (up to 6 missing teeth)
TICML	Titanium AccuraMesh Large (7 or more missing teeth)
PCMS	PEEK AccuraMesh Standard (up to 6 missing teeth)
PCML	PEEK AccuraMesh Large (7 or more missing teeth)

* Please contact your radiologist or device manufacturer if you have any questions on DICOM export.

Design And Order Process

1. Data Submission

Fill in the product request form (online or using form ZB1210) and transfer together with CT/CBCT data (DICOM format required) using Zimmer Biomet Dental's upload website: www.zbd-fileupload.com.

2. Design Phase

The AccuraMesh will be designed according to the requirements written on the request form and you will receive an e-mail with a 3D-PDF file for review. Adjustments can be made to the design at any time prior to final approval.

3. Approval


Once the design is finalized, your approval is required to release the mesh for manufacturing.

4. Manufacturing Of The Mesh

Titanium AccuraMeshes are manufactured by Selective Laser Melting (SLM) and PEEK AccuraMeshes by Fused Deposition Modeling (FDM). The final products are ETO sterilized and provided in a blistered, sterile packaging.

5. Shipment

Once the final design has been approved (step 3), the AccuraMesh will be released from manufacturing after approximately two weeks. The expected delivery date will be confirmed to allow you to schedule surgery accordingly.



Clinical Case




Fig. 1 Left posterior mandible defect.




Fig. 2 Designed Titanium AccuraMesh, lateral view.



Fig. 3 Sterile Titanium AccuraMesh.



Fig. 4 Titanium AccuraMesh fixed in place with osteosynthesis screws. Site grafted with 50:50 mix of autogenous bone and xenograft.



Fig. 5 Titanium AccuraMesh covered with a collagen membrane.



Fig. 6 Radiograph taken after closing the surgical site.



Fig. 7 Soft tissue after 6 months healing time before re-entry.



Fig. 8 Revascularized newly formed bone after mesh removal.



Fig. 9 Prosthetic restoration.




Fig. 10 Radiograph 1-year follow-up.

Clinical photographs courtesy of Dr. N. Cruz. Individual results may vary.

References

1. McAllister BS et al. Journal of Periodontology (2007) 78:377-96.
2. Cruz N et al. Materials (2020) 13.
3. Sghaireen MG et al. Diagnostics (2020) 10:406.
4. Titanium AccuraMesh IFU latest revision.
5. PEEK AccuraMesh IFU latest revision.
6. Her S et al. Journal of Oral and Maxillofacial Surgery (2012) 70:803-810.
7. Parthasarathy J. Ann Maxillofac Surg (2014) 4:9-18.
8. El Chaar E et al. Int J Periodontics Restorative Dent (2019) 39:491-500.
9. Birg J., et al. J Dent Oro Surg (2015) 1:1-4.
10. von Arx T et al. Int J Oral Maxillofac Implants (1996) 11:387-94.
11. Rimondini L et al. J Periodontol (1997) 68:556-62.
12. Ponsonnet L et al. Materials Science and Engineering: C (2002) 21:157-165.
13. Technical Information VESTAKEEP® i4 G PEEK.
14. D'Ercole S et al. Journal of Materials Science: Materials in Medicine (2020) 31:84.
15. Ronda M et al. Int J Periodontics Restorative Dent (2011) 31:505-13.
16. Ronda M et al. Int J Periodontics Restorative Dent (2015) 35:795-801.
17. Heller AL et al. Journal of Oral Implantology (2000) 26:91-103.
18. Romanos GE. Journal of Oral Implantology (2010) 36:25-30.



Contact us at 1-800-342-5454 or visit

zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410, USA
Tel.: +1 561 776 6700
Fax: +1 561 776 1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. AccuraMesh products are manufactured by ResDevMed Lda. Portugal. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZBINST0068 REV A 04/21 ©2021 Zimmer Biomet. All rights reserved.

