





Tissue punches set acc. to Dr. Ady Palti

Ref.-no.: 13388.00

consisting of:

- | | | | |
|----------|-------------------------|---|---|
| 13388.03 | ∅ 3 mm |  | ¼ |
| 13388.04 | ∅ 4 mm |  | |
| 13388.05 | ∅ 5 mm |  | |
| 13388.06 | ∅ 6 mm |  | |
| 6668.01 | Rack for tissue punches | | |



For minimal invasive implant insertion, e.g. in case of one-stage-surgery or for uncovering the implants in case of two-stage-surgery after implant healing period.

