

Instrument for application of bone material with pusher acc. to Prof. Dr. Fouad Khoury

The open applicator for bone material with pusher (piston) has the advantage of enabling the user to pick up and apply bone material. The bone material is picked up by the opening which is designed into one side of the instrument end and can be carried and applied by the piston guide by one hand only via finger. The stodger which is designed on the other end of the instrument allows a condensation and adaption of bone material at the designated defective spot.

Ref.No. **14634.00**



1/1

working end with piston

working end with plugger



piston drawn back



Detail view plugger



with bone graft material