


## Tissue punches set acc. to Dr. Ady Palti


Ref.-no.: **13388.00**

(items zijn ook los te bestellen)

consisting of:

**13388.03**    ø 3 mm    

**13388.04**    ø 4 mm    

**13388.05**    ø 5 mm    

**13388.06**    ø 6 mm    

¼



**6668.02**    Rack for trephine drills and/or tissue punches, 35mmL

For minimal invasive implant insertion, e.g. in case of one-stage-surgery or for uncovering the implants in case of two-stage-surgery after implant healing period.

