



## Prosthetic Overview



# TSX<sup>®</sup>, Trabecular Metal<sup>®</sup>, TSV<sup>®</sup>, and Eztetic<sup>®</sup> Implant Systems



# Trabecular Metal, Eztetic, Tapered Screw-Vent, and TSX Implants

	Length (mm)	Trabecular Metal Implants							
		3.5 mmP				4.5 mmP		5.7 mmP	
		3.7 mmD		4.1 mmD		4.7 mmD		6.0 mmD	
Machined Collar (0.5 mm) with Microgrooves	10		TMMB10		TMM4B10		TMMWB10		TMM6B10
	11,5		TMMB11		TMM4B11		TMMWB11		TMM6B11
	13		TMMB13		TMM4B13		TMMWB13		TMM6B13
	16		TMMB16		-		-		-
MTX Textured Collar with Microgrooves	10		TMTB10		TMT4B10		TMTWB10		TMT6B10
	11,5		TMTB11		TMT4B11		TMTWB11		TMT6B11
	13		TMTB13		TMT4B13		TMTWB13		TMT6B13
	16		TMTB16		-		-		-

	Length (mm)	Eztetic Implants		Tapered Screw-Vent Implants				
		2.9 mmP		3.5 mmP		4.5 mmP		5.7 mmP
		3.1 mmD		3.7 mmD	4.1 mmD	4.7 mmD		6.0 mmD
Machined Collar (1 mm) without Microgrooves	8		TSVB8		TSV4B8		TSVWB8	TSV6B8
	10		TSVB10		TSV4B10		TSVWB10	TSV6B10
	11,5		TSVB11		TSV4B11		TSVWB11	TSV6B11
	13		TSVB13		TSV4B13		TSVWB13	TSV6B13
Machined Collar (0.5 mm) with Microgrooves	16		TSVB16		TSV4B16		TSVWB16	TSV6B16
	8		CM318		TSVM4B8		TSVMWB8	TSVM6B8
	10		CM3110		TSVM4B10		TSVMWB10	TSVM6B10
	11,5		CM3111		TSVM4B11		TSVMWB11	TSVM6B11
MTX Textured Collar with Microgrooves	13		CM3113		TSVM4B13		TSVMWB13	TSVM6B13
	16		CM3116		TSVM4B16		TSVMWB16	TSVM6B16
	8		CT318		TSVT4B8		TSVTWB8	TSVT6B8
	10		CT3110		TSVT4B10		TSVTWB10	TSVT6B10
	11,5		CT3111		TSVT4B11		TSVTWB11	TSVT6B11
	13		CT3113		TSVT4B13		TSVTWB13	TSVT6B13
	16		CT3116		TSVT4B16		TSVTWB16	TSVT6B16

Length (mm)	TSX Implants					
	2.9 mmP	3.5 mmP			4.5 mmP	
	3.1 mmD	3.7 mmD	4.1 mmD	4.7 mmD <sup>1</sup>	5.4 mmD	6.0 mmD <sup>1</sup>
8	TSX31B8	TSX37B8	TSX41B8	TSX47B8	TSX54B8	TSX6B8
10	TSX31B10	TSX37B10	TSX41B10	TSX47B10	TSX54B10	TSX6B10
11,5	TSX31B11	TSX37B11	TSX41B11	TSX47B11	TSX54B11	TSX6B11
13	TSX31B13	TSX37B13	TSX41B13	TSX47B13	TSX54B13	TSX6B13
16	TSX31B16	TSX37B16	TSX41B16	TSX47B16	TSX54B16	TSX6B16

<sup>1</sup> Note: The 4.7 mmD implant features a green color-coded 3.5 mmD platform and the 6.0 mmD implant has a purple color-coded 4.5 mmD platform, unlike TSV Implants.



# Prosthetic Components for TSX, Eztetic, Trabecular Metal, and Tapered Screw-Vent Implants

Height (mm)	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants									
	2.9 mmP		3.5 mmP			4.5 mmP				5.7 mmP <sup>1</sup>		
	3.7 mmEP	4.5 mmEP	3.8 mmEP	5.0 mmEP	6.5 mmEP	4.5 mmEP	5.5 mmEP	6.5 mmEP	7.5 mmEP	5.7 mmEP	6.8 mmEP	7.5 mmEP
3	CEEHA333	CEEHA343	TEEHA333	TEEHA353	TEEHA363	TEEHA443	TEEHA453	TEEHA463	TEEHA473	TEEHA553	TEEHA563	TEEHA573
5	CEEHA335	CEEHA345	TEEHA335	TEEHA355	TEEHA365	TEEHA445	TEEHA455	TEEHA465	TEEHA475	TEEHA555	TEEHA565	TEEHA575
7	CEEHA337	CEEHA347	TEEHA337	TEEHA357	TEEHA367	TEEHA447	TEEHA457	TEEHA467	TEEHA477	-	-	-

One-Piece Healing Collars	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants													
	2.9 mmP		3.5 mmP			4.5 mmP				5.7 mmP <sup>1</sup>						
	Dimension (mm)		Dimension (mm)			Dimension (mm)				Dimension (mm)						
Healing Collars 3 mm length shown	3.7 EP	1.5 CH	CHCNP31		3.5 EP	3 CH	HC333		4.5 EP	3 CH	HC443		6.5 EP	3 CH	HC563	
		3 CH	CHCNP33			5 CH	HC335			5 CH	HC445			5 CH	HC565	
	4.5 EP	4.5 CH	CHCNP34		4.5 EP	3 CH	HC343		5.5 EP	3 CH	HC453		4.5 EP	3 CH	HC563	
		1.5 CH	CHCNP41			5 CH	HC345			5 CH	HC455			5 CH	HC565	
	4.5 EP	3 CH	CHCNP43		5.5 EP	7 CH	HC347		6.5 EP	7 CH	HC457		4.5 EP	3 CH	HC563	
		4.5 CH	CHCNP44			3 CH	HC353			3 CH	HC463			5 CH	HC465	

Impression Copings	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants													
	2.9 mmP		3.5 mmP			4.5 mmP				5.7 mmP <sup>1</sup>						
	Dimension (mm)		Dimension (mm)			Dimension (mm)				Dimension (mm)						
Indirect Transfer (For closed-tray impressions - includes screw p/n HLT52)	3.7 EP	1.5 CH	CITNP31		3.5 EP	3 EP	HLT3/3		4.5 EP	4 EP	HLT4/4		6.5 EP	3 EP	HLT5/6	
		3 CH	CITNP33			4 EP	HLT3/4			5 EP	HLT4/5			5 EP	HLT5/6	
		4.5 CH	CITNP34			5 EP	HLT3/5			6 EP	HLT4/6			6 EP	HLT5/6	
	4.5 EP	1.5 CH	CITNP41		4.5 EP	3 EP	DHT3/3		5.5 EP	4 EP	DHT4/4		6.5 EP	3 EP	DHT5/6 + HLTE	
		3 CH	CITNP43			4 EP	DHT3/4			5 EP	DHT4/5			5 EP	DHT5/6 + HLTE	
4.5 EP	4.5 CH	CITNP44		5.5 EP	5 EP	DHT3/5		6.5 EP	6 EP	DHT4/6		6.5 EP	6 EP	DHT5/6 + HLTE		
	1.5 CH	CDTNP41			3 EP	DHT3/3			4 EP	DHT4/4			5 EP	DHT5/6 + HLTE		
4.5 EP	3 CH	CDTNP43		4.5 EP	4 EP	DHT3/4		5.5 EP	5 EP	DHT4/5		6.5 EP	5 EP	DHT5/6 + HLTE		
	4.5 CH	CDTNP44			5 EP	DHT3/5			6 EP	DHT4/6			6 EP	DHT5/6 + HLTE		
Implant Analogs		CIANP		IA3		IA4		IA5								

Provisional Restorations	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants										
	2.9 mmP		3.5 mmP			4.5 mmP				5.7 mmP <sup>1</sup>			
	Dimension (mm)		Dimension (mm)			Dimension (mm)				Dimension (mm)			
Temporary Titanium/ Plastic Abutments (Plastic Abutments include retaining screw p/n PTSC and long processing screw p/n DHTS2, for max. 180 days in the case of Plastic Abutments)	3.5 EP	CTANP31		PT341S	PT344S	PT341A	PT344A	PT451S	PT454S	PT451A	PT454A	PT561S	PT564S
Fixture Mount/Transfer <sup>2</sup> (Incl. replacement retaining screw p/n URS2)				FMT3				FMT4				FMT5	

<sup>1</sup> Not compatible with TSX Restorative Platform

<sup>2</sup> Supplied with implant

# GenTek® Ti-Bases and Restorative Components

Crowns and Bridges	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants					
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>		
<b>Ti-Base Engaging, Single-Unit, Implant-Level</b> incl. Hexabobular Retaining Screw and Gluing Screw Recommended Torque 30 Ncm Post Height 7 mm								
	1.0 mmH	ZFX11ZB-EZ29AS10E	0.3 mmH	ZFX11ZB-TSV35AS03E	0.3 mmH	ZFX11ZB-TSV45AS03E	0.3 mmH	ZFX11ZB-TSV57AS03E
	1.7 mmH	ZFX11ZB-EZ29AS17E	1.3 mmH	ZFX11ZB-TSV35AS13E	1.3 mmH	ZFX11ZB-TSV45AS13E	1.3 mmH	ZFX11ZB-TSV57AS13E
<b>Ti-Base Non-Engaging, Multi-Unit, Implant-Level</b> incl. Hexabobular Retaining Screw and Gluing Screw Recommended Torque 30 Ncm Post Height 7 mm								
	1.0 mmH	ZFX11ZB-EZ29AS10N	0.3 mmH	ZFX11ZB-TSV35AS03N	0.3 mmH	ZFX11ZB-TSV45AS03N	0.3 mmH	ZFX11ZB-TSV57AS03N
	1.7 mmH	ZFX11ZB-EZ29AS17N	1.3 mmH	ZFX11ZB-TSV35AS13N	1.3 mmH	ZFX11ZB-TSV45AS13N	1.3 mmH	ZFX11ZB-TSV57AS13N
<b>Hexabobular Retaining Screw for Ti-Base</b>								
	ZFX09-ZB-EZ-HLRS	ZFX09-ZB-TSV-HLRS	ZFX09-ZB-TSV-HLRS	ZFX09-ZB-TSV-HLRS				
<b>Ti-Base Milling Block Compatible Engaging, Single-Unit, Implant-Level</b> incl. Retaining and Gluing Screws Recommended Torque 30 Ncm Height 4.7 mm								
	ZFX11-ZB-EZ-2947-ES	ZFX11-ZB-TSV-3547-EL	ZFX11-ZB-TSV-4547-EL	ZFX11-ZB-TSV-5747-EL				
<b>Ti-Base Milling Block Compatible Non-Engaging, Multi-Unit, Implant-Level</b> incl. Retaining and Gluing Screws Recommended Torque 30 Ncm Height 4.7 mm								
	ZFX11-ZB-EZ-2947-NES	ZFX11-ZB-TSV-3547-NEL	ZFX11-ZB-TSV-4547-NEL	ZFX11-ZB-TSV-5747-NEL				
<b>Retaining Screw for Ti-Base Milling Block Compatible</b>								
	ZFX09-ZB-EZ-RS	ZFX09-ZB-TSV-RS	ZFX09-ZB-TSV-RS	ZFX09-ZB-TSV-RS				
<b>Ti-Base Non-Engaging, Multi-unit, Abutment-Level</b> incl. Retaining and Gluing Screws Recommended Torque 20 Ncm Height 3.9 mm								
			ZFX11-ZB-TA-45					
<b>Retaining Screw for Multi-unit Ti-Base, Abutment-Level</b>								
			ZFX09-ZB-TA-RS					

Pre-Abutment Blanks	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
<b>Pre-Abutment Blank (Titanium)</b> incl. Retaining Screw, 30 Ncm Recommended Torque 30 Ncm						
	ZFX18-ZB-EZ-29	ZFX18-ZB-TSV-35	ZFX18-ZB-TSV-45	ZFX18-ZB-TSV-57		
<b>Retaining Screw for Ti-Base</b>						
	ZFX09-ZB-EZ-RS	ZFX09-ZB-TSV-RS	ZFX09-ZB-TSV-RS	ZFX09-ZB-TSV-RS		

Scanbodies	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
<b>Desktop Scanbody, Implant-Level</b> incl. Screw						
	ZFX05D-ZB-EZ-29	ZFX05D-ZB-TSV-35	ZFX05D-ZB-TSV-45	ZFX05D-ZB-TSV-57		
<b>Intraoral Scanbodies, Implant-Level</b> incl. Screws, Double pack Available in heights of 4.0 mm and 7.0 mm						
	ZFX05I-ZB-EZ-29	ZFX05I-ZB-TSV-35	ZFX05I-ZB-TSV-45	ZFX05I-ZB-TSV-57		
<b>Desktop Scanbody, Abutment Level</b> incl. Screw						
			ZFX05D-ZB-TA-45			
<b>Intraoral Scanbody, Abutment-Level</b> incl. Screws, Double pack Available in height of 4.0 mm only						
			ZFX05I-ZB-TA-45			

<sup>1</sup> Not compatible with TSX Restorative Platform

Digital Scan Analogs	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
Digital Scan Analog, Implant-Level	 ZFX15-ZB-EZ-29		 ZFX15-ZB-TSV-35		 ZFX15-ZB-TSV-45	 ZFX15-ZB-TSV-57
Placement Tool for Digital Scan Analogs	 ZFX02-M16		 ZFX02-M18		 ZFX02-M18	 ZFX02-M18
Digital Scan Analogs, Abutment Level	-		-		 ZFX15-ZB-TA-45	-
Placement Tool for Digital Scan Analog Abutment Level	-		-		 ZFX02-M18	-

## Cement-Retained Restorations

Straight	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
Hex-Lock® Contour Abutments (Incl. screw p/n MHLAS)	Dimension (mm)  3.7 EP 1.5 CH <b>CANP31S</b> 3 CH <b>CANP33S</b> 4.5 CH <b>CANP34S</b> 4.5 EP 1.5 CH <b>CANP41S</b> 3 CH <b>CANP43S</b> 4.5 CH <b>CANP44S</b>		Dimension (mm)  4.5 EP 1 CH <b>ZOA341S</b> 2 CH <b>ZOA342S</b> 3 CH <b>ZOA343S</b>		Dimension (mm)  4.5 EP 1 CH <b>ZOA441S</b> <sup>2</sup> 2 CH <b>ZOA442S</b> <sup>2</sup> 3 CH <b>ZOA443S</b> <sup>2</sup> 5.5 EP 1 CH <b>ZOA451S</b> 2 CH <b>ZOA452S</b> 3 CH <b>ZOA453S</b>	Dimension (mm)  6.5 EP 1 CH <b>ZOA561S</b> 2 CH <b>ZOA562S</b> 3 CH <b>ZOA563S</b>
Contour Restorative Kits (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)	-		4.5 EP <b>ZCRKIT4S</b> 		5.5 EP <b>ZCRKIT5S</b> 	6.5 EP <b>ZCRKIT6S</b> 
Hex-Lock Short Abutments (Incl. screw p/n MHLAS)	-		 4.5 EP 1 CH <b>SA341</b> 2 CH <b>SA342</b>		 5.5 EP 1 CH <b>SA451</b> 2 CH <b>SA452</b>	 6.5 EP 1 CH <b>SA561</b> 2 CH <b>SA562</b>
Short Restorative Kits (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)	-		4.5 EP <b>SARKIT4</b> 		5.5 EP <b>SARKIT5</b> 	6.5 EP <b>SARKIT6</b> 
Hex-Lock Abutments (Incl. screw p/n MHLAS)	-		 3.5 EP <b>HLA3/3</b> 4.5 EP <b>HLA3/4</b> 5.5 EP <b>HLA3/5</b>		 4.5 EP <b>HLA4/4</b> 5.5 EP <b>HLA4/5</b> 6.5 EP <b>HLA4/6</b>	 6.5 EP <b>HLA5/6</b>

Angled	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
Hex-Lock Contour Abutments, 17° (Incl. screw p/n AH20S)	Dimension (mm)  3.7 EP 1.5 CH <b>CANP31A</b> 3 CH <b>CANP33A</b> 4.5 EP 1.5 CH <b>CANP41A</b> 3 CH <b>CANP43A</b>		Dimension (mm)  4.5 EP 1 CH <b>ZOA341A</b> 2 CH <b>ZOA342A</b>		Dimension (mm)  4.5 EP 1 CH <b>ZOA441A</b> <sup>2</sup> 2 CH <b>ZOA442A</b> <sup>2</sup> 5.5 EP 1 CH <b>ZOA451A</b> 2 CH <b>ZOA452A</b>	Dimension (mm)  6.5 EP 1 CH <b>ZOA561A</b> 2 CH <b>ZOA562A</b>
Contour Restorative Kits, 17° (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)	-		4.5 EP <b>ZCRKIT4S</b> 		5.5 EP <b>ZCRKIT5S</b> 	6.5 EP <b>ZCRKIT6S</b> 
Angled Abutment, 20°, for 6 Positions (Incl. screw p/n AH20S)	 3.5 EP <b>C20A3</b>		 4.5 EP <b>AH20/4</b>		 5.5 EP <b>AH20W/5</b>	 6.5 EP <b>A5H20/6</b>

<sup>1</sup> Not compatible with TSX Restorative Platform

<sup>2</sup> Structures for Ø 4.5 mm x Ø 4.5 mm profile are listed but not illustrated. Note: For these items, please use the Contour Restorative Kit ZCRKIT4S (straight) or ZCRKIT4A (angled).

## Custom Restorations, Engaging and Non-Engaging

Custom Restorations	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
“Cast-To” Gold Abutment, Engaging (Incl. screw p/n MHLAS)	-		HLA3G	HLA4G	HLA5G	
“Cast-To” Gold Abutment, Non-Engaging (Incl. screw p/n MHLAS)	CEANP31		NEA3G	NEA4G	-	

## Replacement Retaining Screws

Screws	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants			
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>
<b>Abutment Screw</b> For Straight Abutments	CUAS		MHLAS		MHLAS	MHLAS
<b>Abutment Screw</b> For Angled Abutments	CUASA 1.5 mmCH	CUAS 3 mmCH	AH20S		AH20S	AH20S
<b>Replacement Retaining Screw</b> Direct Transfer	CASLC		DHTS		DHTS	DHTS
<b>Replacement Retaining Screw</b> Plastic Temporary Abutments (Long)	-		DHTS2		DHTS2	DHTS2
<b>Extension Screw</b> 3.0 mm, Indirect Transfer	-		HLTE		HLTE	HLTE
<b>Replacement Retaining Screw</b> Indirect Transfer, Hex-Lock Abutments (Long)	-		HLTS2		HLTS2	HLTS2
<b>Waxing Screw</b> “Cast-To” Gold Abutments, Engaging	-		MTWSD		MTWSD	MTWSD
<b>Waxing Screw</b> Tapered Abutment	-		SCWS		SCWS	SCWS
<b>Replacement Screw For Copings</b>	-		SCTS		SCTS	SCTS

## Screw-Retained Restorations

Universal Multi-Unit Abutments	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants				
	2.9 mmP		3.5 mmP		4.5 mmP	5.7 mmP <sup>1</sup>	
<b>Universal Straight Abutment</b> Non-engaging, does not engage in the internal hex connection. Not for use in single-unit restorations.	-		1.0 mmH  UMUA-TSV351	2.0 mmH  UMUA-TSV352	3.0 mmH  UMUA-TSV353	4.0 mmH  UMUA-TSV354	5.0 mmH  UMUA-TSV355
<b>Universal Angled Abutment, 17°</b> Engaging in implant, not for use in single-unit restorations (incl. preassembled retaining screw)	-		2.5 mmH  UMUA-TSV3525-17	3.5 mmH  UMUA-TSV3535-17	4.5 mmH  UMUA-TSV4525-17	5.5 mmH  UMUA-TSV4535-17	-
<b>Universal Angled Abutment, 30°</b> Engaging in implant, not for use in single-unit restorations (incl. preassembled retaining screw)	-		3.5 mmH  UMUA-TSV3535-30	4.5 mmH  UMUA-TSV3545-30	5.5 mmH  UMUA-TSV4535-30	6.5 mmH  UMUA-TSV4545-30	-

### Universal Multi-Unit Abutment Components



<sup>1</sup> Not compatible with TSX Restorative Platform

# Overdenture Restorations

OverdenSURE® Abutments	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants					
	2.9 mmP		3.5 mmP		4.5 mmP		5.7 mmP <sup>1</sup>	
<b>OverdenSURE</b> Removable Attachment System (Retention Kit ODS-RIKIT2001 included. Available in all active geographies except the U.S.)			1.0 mmH	ODS-TSV301	1.0 mmH	ODS-TSV401	1.0 mmH	ODS-TSV501
			2.0 mmH	ODS-TSV302	2.0 mmH	ODS-TSV402	2.0 mmH	ODS-TSV502
			3.0 mmH	ODS-TSV303	3.0 mmH	ODS-TSV403	3.0 mmH	ODS-TSV503
			4.0 mmH	ODS-TSV304	4.0 mmH	ODS-TSV404	4.0 mmH	ODS-TSV504
			5.0 mmH	ODS-TSV305	5.0 mmH	ODS-TSV405	5.0 mmH	ODS-TSV505
			6.0 mmH	ODS-TSV305	6.0 mmH	ODS-TSV406	6.0 mmH	ODS-TSV506

OverdenSURE Components														
Transfer Cap		QTY	Abutment Analogue		QTY	Bar Attachment <sup>2</sup>		QTY	Blockout Spacer		QTY	Retention Housing Cap		QTY
	ODS-IC01	1		ODS-AA01	1		ODS-BARATTO1	1		ODS-BS010	20		ODS-RH01	1
	ODS-IC04	4		ODS-AA04	4		ODS-BARATTO2	2					ODS-RH04	4
	ODS-IC10	10		ODS-AA10	10		ODS-BARATTO4	4					ODS-RH10	10
Retention Kits (Set)				QTY	Retention Inserts				QTY					
	20°	ODS-RIKIT2001	1		ODS-RIB2004	4		ODS-RIP2004	4		ODS-RIR4004	4		
	20°	ODS-RIKIT2002	2		ODS-RIP2004	4		ODS-RIR4004	4		ODS-RIY4004	4		
	40°	ODS-RIKIT4001	1		ODS-RIC2004	4		ODS-RIG4004	4					
	40°	ODS-RIKIT4002	2											

Ball Abutments	TSX and Eztetic Implants		Trabecular Metal, Tapered Screw-Vent and TSX Implants					
	2.9 mmP		3.5 mmP		4.5 mmP		5.7 mmP <sup>1</sup>	
<b>Ball Abutments</b> with Metal Housing and Nylon Liners	2 mmCH	CBANP2	2 mmCH	BAC2	2 mmCH	BACW2	2 mmCH	BA5C2
	4 mmCH	CBANP4	4 mmCH	BAC4	4 mmCH	BACW4	4 mmCH	BA5C4
	6 mmCH	CBANP6	6 mmCH	BAC6	6 mmCH	BACW6	6 mmCH	

Transfer Components for Ball Abutments					
	<b>BAT</b> Ball Abutment Transfer Cap		<b>CAT</b> Transfer Cap		<b>BAR</b> Ball Abutment Analogue
	<b>CAH</b> Metal Females		<b>CAN</b> Nylon Liner		<b>CAN-G</b> Nylon Liner for Fixed Retention

Prosthetic Tools and Torque Wrenches					
	<b>TWR</b> Torque Wrench 10-35 Ncm		<b>TW1.25</b> Hex Driver Tip for Torque Wrench Ø 1.25 mm, for TWR, 17 mm long		<b>TW1.25L</b> Hex Driver Tip for Torque Wrench Ø 1.25 mm, long, for TWR, 22 mm long
	<b>HXGR1.25</b> 22 mm long		<b>Hx1.25D</b> 23 mm long		<b>TLRT2</b> Removal Tool TSX
	<b>ZTIRW</b> High Torque Indicating Ratchet Wrench 40-90 Ncm		<b>ABTH</b> Abutment Holder		<b>ODS-DRVR</b> Abutment Driver
	<b>OVS-IRTOOL</b> Removal and Insertion Tool		<b>HXLGR1.25</b> 30 mm long		<b>HXL1.25D</b> 30 mm long
	<b>UMUA-D</b> Multi-Unit Driver		<b>ISO4X4A</b> ISO Latch, 4x4 Driver		<b>ASC-SD24</b> Angulated Screw Channel Screwdriver, 24 mm, long
	<b>ASC-SD30</b> Angulated Screw Channel Screwdriver, 30 mm, long				



schacht 12,5mm, totaal 24mm, TASCD24

schacht 15mm, totaal 30mm, TASCD30

<sup>1</sup> Not compatible with TSX Restorative Platform

<sup>2</sup> Does not include Retention Kit sets



[www.implacom.nl](http://www.implacom.nl) | +31 (0)577 46 1927

[For more information visit ZimVie.com](http://ZimVie.com)

**ZimVie**

4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272  
[dentalCS@ZimVie.com](mailto:dentalCS@ZimVie.com)

**Biomet 3i (UK & Ireland) Ltd**

Reading Business Centre  
Suite 807, 8<sup>th</sup> Floor Fountain House  
2 Queens Walk, Reading, Berks  
RG1 7QF, United Kingdom  
Phone (UK): 0800 652 1233  
Phone (Ireland): 1 800 552 752  
[ZV.UKorders@ZimVie.com](mailto:ZV.UKorders@ZimVie.com)



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. OverdenSURE and Universal Multi-unit abutments, components and tools are manufactured by Terrats Medical S.L. Ratchet is manufactured by Josef Ganter. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie.



ZV0796 REV B 12/23 ©2023 ZimVie. All rights reserved.