




Zfx™ GenTek™

Genuine Restorations and
Open Digital Workflows

Product Catalog



GenTek™ Open Connectivity

For even more flexibility.




Genuine

GenTek™ Pre-Abutment Blanks are compatible with Medentika® holders. (CAM update required!)



Compatible

Compatibility of the GenTek™ Ti-Base to Sirona® CEREC® Blocks



exocad

 **dental wings**

3shape 



Free

Libraries for the GenTek restorative components are available for Zfx™ CAD, exocad, 3Shape and Dentalwings CAD software. Free download at zfx-dental.com

Coming soon in 2019:

- Ti-Bases and ancillary components for Certain®, External Hex® and Eztec® connections
- Ti-Bases for full-arch restorations
- Pre-milled abutment blanks for all 4 connections



SCAN & INFORM
More at
www.zfx-dental.com

 **ZIMMER BIOMET**
A COMPANY OF

Table of Contents

Zimmer Biomet Dental Tapered Screw-Vent® and Trabecular metal connection 4

Zimmer Biomet Dental Eztetic® Connection 6

Placement Tool Compatibility Guide



Placement Tool	Laboratory Analog
M 16	Certain & Eztetic Lab Analogs
M 18	TSV, TM & Tapered Abutment Lab Analogs

Legend

Ø Restorative Platform (mm)
Color code for Eztetic® and Tapered Screw-Vent®

- 2.9mmD
- 3.5mmD
- 4.5mmD
- 5.7mmD

Zimmer Biomet Dental Tapered Screw-Vent® and Trabecular Metal™ Connection

Single-Unit Restorative Components

Ti-Base Engaging¹ *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-EL
	4.5	ZFX11-ZB-TSV-4547-EL
	5.7	ZFX11-ZB-TSV-5747-EL

Retaining Screw



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

Titanium Pre-milled Abutment Blank² *incl. Retaining Screw*



Restorative Platform (mm)	Part No.
3.5	ZFX18-ZB-TSV-35
4.5	ZFX18-ZB-TSV-45
5.7	ZFX18-ZB-TSV-57

Retaining Screw



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

Multi-Unit Restorative Components

Ti-Base Non-Engaging¹ *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-NEL
	4.5	ZFX11-ZB-TSV-4547-NEL
	5.7	ZFX11-ZB-TSV-5747-NEL

Retaining Screw



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

Components Accessories & Tools

Desktop Scanbody incl. Screw



Torque Wrench for Desktop Scanbody



Digital Scan Analog⁴

Placement Tool for Digital Scan Analogs



Restorative Platform (mm)	Part No.	Hand-tighten ³	Part No.	Restorative Platform (mm)	Part No.	Part No.
3.5	ZFX05D-ZB-TSV-35		ZFX02002067	3.5	ZFX15-ZB-TSV-35	ZFX02-M18
4.5	ZFX05D-ZB-TSV-45			4.5	ZFX15-ZB-TSV-45	
5.7	ZFX05D-ZB-TSV-57			5.7	ZFX15-ZB-TSV-57	

Intraoral Scanbodies incl. Screws



Latch-Lock Hex Driver⁵ for Abutment Screws



Height (mm)	Restorative Platform (mm)	Part No.	Description	Part No.
4.0 & 7.0	3.5	ZFX05I-ZB-TSV-35	Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mm D, 23 mm L (May require Reduction Handpiece)	HX1.25D
	4.5	ZFX05I-ZB-TSV-45	Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mm D, 26 mm L (May require Reduction Handpiece)	HXL1.25D
	5.7	ZFX05I-ZB-TSV-57		

Hex Driver for Abutment Screws & Intraoral Scanbodies



Removal Tool



Torque Wrench Tool for Abutment Screws & Intraoral Scanbodies



Description	Part No.	Description	Part No.	Description	Part No.
Hex Driver, Short, with GemLock(R) Retention, 1.25 mm D, 22 mm L	HXGR1.25	Removal Tool for TVS/TM Abutments	TLRT2	Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Hex Driver, Long, with GemLock(R) Retention, 1.25 mm D, 30 mm L	HXLGR1.25	Removal Screw for TSV/TM Abutments	HLRS	Torque Wrench Driver, Short	TW1.25 ⁵
				Torque Wrench Driver, Long	TW1.25L ^{5/6}

¹Groove indicating 4 mm post height

²Coming soon

³Minimum of 5 Ncm; Maximum of 15 Ncm.

Refer to Instruction Manual for more information.

⁴Compatible with existing TSV/ TM Impression Copings

⁵Compatible with Tibases

⁶Compatible with IOS Scanbodies

Zimmer Biomet Dental Eztetic® Connection

Single-Unit Restorative Components

Ti-Base Engaging¹ *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-ES

Retaining Screw



Part No.	
ZFX09-ZB-EZ-RS	30 Ncm

Titanium Pre-milled Abutment Blank *incl. Retaining Screw*



Restorative Platform (mm)	Part No.
2.9	ZFX18-ZB-EZ-29

Retaining Screw



Part No.	
ZFX09-ZB-EZ-RS	30 Ncm

Multi-Unit Restorative Components

Ti-Base Non-Engaging¹ *incl. Retaining and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-NES

Retaining Screw



Part No.	
ZFX09-ZB-EZ-RS	30 Ncm

Components Accessories & Tools

Desktop Scanbody <i>incl. Screw</i>			Torque Wrench <i>for Desktop Scanbody</i>		Digital Scan Analog ³	Placement Tool <i>for Digital Scan Analogs</i>
						
Restorative Platform (mm)	Part No.		Part No.		Restorative Platform (mm)	Part No.
2.9	ZFX05D-ZB-EZ-29	Hand-tighten ²	ZFX02002067		2.9	ZFX15-ZB-EZ-29
						ZFX02-M16

Intraoral Scanbodies <i>incl. Screws</i>				Latch-Lock Hex Driver ⁴ <i>for Abutment Screws</i>	
					
Height (mm)	Restorative Platform (mm)	Part No.		Description	Part No.
4.0 & 7.0	2.9	ZFX05I-ZB-EZ-29	Hand-tighten ²	Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mm D, 23 mm L (May require Reduction Handpiece)	HX1.25D
				Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mm D, 26 mm L (May require Reduction Handpiece)	HXL1.25D

Hex Driver <i>for Abutment Screws & Intraoral Scanbodies</i>		Removal Tool		Torque Wrench Tool <i>for Abutment Screws & Intraoral Scanbodies</i>	
					
Description	Part No.	Description	Part No.	Description	Part No.
Hex Driver, Short, with GemLock(R) Retention, 1.25 mm D, 22 mm L	HXGR1.25	for Eztec Abutments	CLRT2	Torque Wrench, Restorative (adjustable torque range 10–35 Ncm)	TWR
Hex Driver, Long, with GemLock(R) Retention, 1.25 mm D, 30 mm L	HXLGR1.25			Torque Wrench Driver, Short	TW1.25 ⁴
				Torque Wrench Driver, Long	TW1.25L ^{4 5}

¹ Groove indicating 4 mm post height

² Eztec connection Ti-Base require the use of the CEREC blocks for Certain 3.4 mm Connection

³ Minimum of 5 Ncm; Maximum of 15 Ncm. Refer to Instruction Manual for more information.

⁴ Compatible with existing Eztec Impression Copings

⁵ Compatible with Tibases

⁶ Compatible with IOS Scanbodies



To order genuine
GenTek™ Components,
call: +34 93 470 55 00

Zfx GmbH

Kopernikusstraße 15
85221 Dachau, Germany
Tel. +49 81 31 33 244 - 0
Fax +49 81 31 33 244 - 10
www.zfx-dental.com

All content herein is protected by copyright, trademarks or other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated. ©2019 Zimmer Biomet. CEREC and Sirona are trademarks of Dentsply Sirona or its affiliates.

