

BTI EXTRACTORS COMPATIBILITY CHART

IMPLANT BRAND	MODEL	✓	Nº
ANKYLOS	---	✓	3
ASTRA TECH	SMALL (Ø3,5/4)	✓	0
	LARGE (Ø4,5/5)	✓	3
	X-SMALL (Ø3)	✓	1A
AVINENT	HE Ø4	✓	1
BEGO SEMADOS	S Line Ø 3,25	✓	1 - 2
	S Line Ø3.75	✓	3
	S Line Ø4.1	✓	1
BICON	Narrow Implants Ø 3,5	✓	1
	Narrow Implants Ø 4	✓	1
	Short implants Ø 4,5	✓	3
	Short implants Ø 5	✓	3
	Short implants Ø 6	✓	3
BIOHORIZONS	---	✓	1
BIOMEDICS IMPLANTS	---	✓	1
BIOMET 3i	Certain*	✓	1A
BREDENT DENTAL	MINI 2 SKY IMPLANT SYSTEM	✓	0
	BLUE SKY IMPLANT SYSTEM	✓	1 - 2
BTI	INTERNA	✓	3
	INTERNA ESTRECHA	✓	2
	EXTERNA	✓	1
	TINY	✓	1
CAMLOG	SCREW LINE Ø4,3	✓	2
	CAMLOG ROOT LINE Ø 5	✓	4
	CAMLOG SCREW LINE Ø 5	✓	4
	CONELOG SCREW LINE Ø 5	✓	4
DENTAURUM IMPLANTS	TioLogic	✓	0 - 1
	TioLox Ø4,5	✓	0
DENTSPLY FRIADENT	XIVE* S D3	✓	1 - 2
	XIVE* S D3,8	✓	2
	XIVE* S D5,5	✓	3 - 4
	ANKYLOS* C/X A Ø 3,5	✓	1
	ANKYLOS* C/X B Ø 4,5	✓	1
	ANKYLOS* C/X C Ø 5,5	✓	1
	ANKYLOS* C/X D Ø 7	✓	1
	FRIALIT* D3,8	✓	2
	FRIALIT* D4,5	✓	2
	FRIALIT* D5,5	✓	2
	FRIALIT* D6,5	✓	2
ECKERMANN	Implante Roscado DH Ø3,5	✓	1
IMPLANT DIRECT LLC	Screw Direct Ø 3,7	✓	1
	Screw Direct Ø 4,7	✓	1
	Dual Grip Ø4,3	✓	2
	Screw Plant Ø3,7	✓	1
	Screw Plant Ø4,7	✓	1
IMZ	---	✓	1
KLOCKNER	S3M 3112	✓	0
	Cone Essential	✓	3
	Cone 1.5 Essential	✓	3

✓ OPTIMUM COMPATIBILITY TESTED BY BTI

✓ OPTIMUM COMPATIBILITY IN EXPERIMENTAL CLINICAL TRIALS

Nº 0 - EXIME

Nº 1 - EXIM

Max
70 Ncm of torque



Nº 2 - INEXIML

Nº 3 - INEXIM



Nº 4 - EXIVA



Nº 1A - EXIML-A

Nº 1B - EXIML-B



Ed.3/Rev.1 - ENG

IMPLACOM

ZimVie

BTI EXTRACTORS COMPATIBILITY CHART

IMPLANT BRAND	MODEL	✓	Nº
LIFECORE	RENOVA™ SCI	✓	1
	RENOVA™ RDI	✓	1
MEDENTICA	Microcone Implantat Ø4,5/3,5	✓	0
	Microcone Implantat Ø4	✓	0
	M-Implantat Ø3	✓	0
		✓	0
MEGAGEN	Any Ridge Ø4	✓	3
MIS	SEVEN 4,2	✓	2-3
MOZOGRAU	MG-OSSEOUS std Ø 3,4	✓	1
	MG-OSSEOUS std Ø 3,75	✓	1
	MG-OSSEOUS std Ø 4,25	✓	1
NEODENT	CM	✓	1
	GT	✓	1
	Titamax Poros Ø 3,3	✓	1
	Titamax Alvim Ø 3,5	✓	1
	Titamax TI Ø 3,75 – 4,0	✓	1
	Alvim TI Ø 4,3	✓	1
		✓	1
NOBEL BIOCARE	Replace Ø3,5	✓	2
	Replace Ø4,3	✓	2-3
	Replace Ø5	✓	3-4
	Branemark System NP	✓	1
	Branemark System RP	✓	1
	NobelDirect Groovy RP	✓	2
	NobelDirect 3,0	✓	1
		✓	1
PHIBO	TSA ADVANCE S3	✓	0-1
	TSA ADVANCE S4	✓	1
	TSA ADVANCE S5	✓	1
	TSA S3	✓	1
	TSA S4	✓	1
	TSH S3	✓	1
	TSH S4	✓	1
	BNT S2	✓	1*
		✓	1
		✓	1
SOUTHERN IMPLANTS	IA Ø4,3	✓	2
	BAT	✓	1
STRAUMANN	Bone Level (NC)	✓	1
	Regular Neck (RN)	✓	3
TITANIMPLANT	Serie Avanza Ø3,3	✓	1
TRS	---	✓	1**
ZIMMER	Tapered Screw Vent	✓	1A

Note: THOMMEN implantaten: de maatvoering van Zimmer is ook te gebruiken voor het Thommen systeem.

(*) This 3.3 diameter Implant can not resist the torque applied by the LL7200 wrench so the LL400 surgical wrench should be used in this case.
(**) 1 cut by the second thread.

This document shows confidential information and it can only be in BTI employees hands. It is absolutely forbidden to deliver, distribute, publish or any other misuse of it, if so it will be considered a serious breach of the non disclosure agreement and BTI will be entitled to take any legal actions in its hands to defend its rights.

✓ OPTIMUM COMPATIBILITY TESTED BY BTI
✓ OPTIMUM COMPATIBILITY IN EXPERIMENTAL CLINICAL TRIALS

Nº 0 - EXIME

Máx.
70 Ncm of torque



Nº 1 - EXIM



Nº 2 - INEXIML



Nº 3 - INEXIM



Nº 4 - EXIMA



Nº 1A - EXIML-A

Nº 1B - EXIML-B

